

Hydrogen Peroxide, 30%
BAKER
(Stabilized)
Formerly C.P.



Material No.: 2189-01
Batch No.: 0000014468
Manufactured Date: 2012/08/10
Expiration Date: 2014/02/08

Certificate of Analysis

Test	Specification	Result
Assay (H ₂ O ₂)	29.0 - 32.0 %	31.5
Residue after Evaporation	<= 0.0300 %	< 0.0005
Nitrogen Compounds (as N)	<= 0.001 %	< 0.001
Phosphate (PO ₄)	<= 0.02 %	< 0.01
Sulfate (SO ₄)	<= 0.001 %	< 0.001
Chloride (Cl)	<= 5 ppm	< 1
Heavy Metals (as Pb)	<= 1 ppm	< 1
Trace Impurities - Iron (Fe)	<= 0.500 ppm	< 0.100

For Laboratory, Research, or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC

ISO

Phillipsburg, NJ 9001.2008, 14001.2004
Paris, KY 9001.2008
Mexico city, Mexico 9001.2008
Deventer, The Netherlands 9001.2008, 14001.2004
Selangor, Malaysia 9001.2008
Panoli, India 9001.2008
Gliwice, Poland 9001.2008, 17025.2005

Richard M Siberski
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.
3477 Corporate Parkway, Suite #200, Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610