

Hydrogen Peroxide, 30%
BAKER
(Stabilized)
Formerly C.P.



Material No.: 2189-01
Batch No.: 0000002047
Manufactured Date: 2012/04/06
Expiration Date: 2013/10/05

Certificate of Analysis

Test	Specification	Result
Assay (H ₂ O ₂)	29.0 - 32.0 %	31.1
Residue after Evaporation	<= 0.03 %	0.00
Nitrogen Compounds (as N)	<= 0.001 %	<0.001
Phosphate (PO ₄)	<= 0.02 %	<0.00
Sulfate (SO ₄)	<= 0.001 %	<0.000
Chloride (Cl)	<= 5 ppm	<0
Heavy Metals (as Pb)	<= 1 ppm	<0
Trace Impurities - Iron (Fe)	<= 0.5 ppm	<0.0

For Laboratory, Research, or Manufacturing Use

Storage Conditions:

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC

ISO

Phillipsburg, NJ 9001.2008, 14001.2004
Paris, KY 9001.2008
Mexico city, Mexico 9001.2008
Deventer, The Netherlands 9001.2008, 14001.2004
Selangor, Malaysia 9001.2008
Panoli, India 9001.2008
Gliwice, Poland 9001.2008, 17025.2005

Richard M Siberski
Global Director of Quality Assurance