

2-Propanol
HPLC
For use in Liquid Chromatography (HPLC & UHPLC) &
Spectrophotometry



Material No.: 9095-33
Batch No.: 0000001440
Manufactured Date: 2012/04/20

Certificate of Analysis

Test	Specification	Result
Ultraviolet Absorbance (1.00-cm cell v	<= 0.01	< 0.00
Ultraviolet Absorbance (1.00-cm cell v	<= 0.01	0.00
Ultraviolet Absorbance (1.00-cm cell v	<= 0.020	0.006
Ultraviolet Absorbance (1.00-cm cell v	<= 0.16	0.06
Ultraviolet Absorbance (1.00-cm cell v	<= 205	203
Assay (CH ₃ CHOHCH ₃)	>= 99.7 %	100.0
Residue after Evaporation	<= 2 ppm	< 1
Water (by KF, coulometric)	<= 0.03 %	0.01

For Laboratory, Research, or Manufacturing Use

Storage Conditions:

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC

ISO

Phillipsburg, NJ 9001.2008, 14001.2004
Paris, KY 9001.2008
Mexico city, Mexico 9001.2008
Deventer, The Netherlands 9001.2008, 14001.2004
Selangor, Malaysia 9001.2008
Panoli, India 9001.2008
Gliwice, Poland 9001.2008, 17025.2005

Richard M Siberski
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.
3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610