

2-Propanol  
HPLC  
For use in Liquid Chromatography (HPLC & UHPLC) &  
Spectrophotometry



Material No.: 9095-33  
Batch No.: 0000003501  
Manufactured Date: 2012/05/15

## Certificate of Analysis

Test	Specification	Result
Ultraviolet Absorbance (1.00-cm cell v	$\leq 0.01$	< 0.00
Ultraviolet Absorbance (1.00-cm cell v	$\leq 0.01$	0.00
Ultraviolet Absorbance (1.00-cm cell v	$\leq 0.020$	0.006
Ultraviolet Absorbance (1.00-cm cell v	$\leq 0.16$	0.08
Ultraviolet Absorbance (1.00-cm cell v	$\leq 205$	203
Assay (CH <sub>3</sub> CHOHCH <sub>3</sub> )	$\geq 99.7\%$	99.9
Residue after Evaporation	$\leq 2.0$ ppm	< 1.0
Water (by KF, coulometric)	$\leq 0.03\%$	0.02

For Laboratory, Research, or Manufacturing Use

Storage Conditions:

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC

ISO

Phillipsburg, NJ 9001.2008, 14001.2004  
Paris, KY 9001.2008  
Mexico city, Mexico 9001.2008  
Deventer, The Netherlands 9001.2008, 14001.2004  
Selangor, Malaysia 9001.2008  
Panoli, India 9001.2008  
Gliwice, Poland 9001.2008, 17025.2005

Richard M Siberski  
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.  
3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610