

2-Propanol
HPLC
For use in Liquid Chromatography (HPLC & UHPLC) &
Spectrophotometry



Material No.: 9095-33
Batch No.: 0000008979
Manufactured Date: 2012/07/04
Retest Date: 2017/07/03

Certificate of Analysis

Test	Specification	Result
Ultraviolet Absorbance (1.00-cm cell vs. water) - 350 nm	≤ 0.01	<0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) - 280 nm	≤ 0.01	<0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) - 254 nm	≤ 0.020	0.013
Ultraviolet Absorbance (1.00-cm cell vs. water) - 225 nm	≤ 0.16	0.10
Ultraviolet Absorbance (1.00-cm cell vs. water) - UV Cut-off, nm	≤ 205	203
Assay (CH ₃ CHOHCH ₃)	$\geq 99.7\%$	99.8
Residue after Evaporation	≤ 2.0 ppm	<1.0
Water (by KF, coulometric)	$\leq 0.03\%$	0.01

For Laboratory, Research, or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC

ISO

Phillipsburg, NJ 9001.2008, 14001.2004
Paris, KY 9001.2008
Mexico city, Mexico 9001.2008
Deventer, The Netherlands 9001.2008, 14001.2004
Selangor, Malaysia 9001.2008
Panoli, India 9001.2008
Gliwice, Poland 9001.2008, 17025.2005

Richard M Siberski
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.
3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610