

1-Octanesulfonic Acid, Sodium Salt  
HPLC  
For Ion-Pair Chromatography of Basic Compounds



Material No.: 2818-05  
Batch No.: 0000049796  
Manufactured Date: 2013/05/01  
Retest Date: 2020/04/29

## Certificate of Analysis

Test	Specification	Result
Assay (as C <sub>8</sub> H <sub>17</sub> SO <sub>3</sub> Na) (by acidimetry)	>= 98.0 %	99.8
UV Absorbance of 0.25 M Solution - 250 nm	<= 0.050	0.008
UV Absorbance of 0.25 M Solution - 240 nm	<= 0.050	0.012
UV Absorbance of 0.25 M Solution - 230 nm	<= 0.050	0.020
UV Absorbance of 0.25 M Solution - 220 nm	<= 0.060	0.035
UV Absorbance of 0.25 M Solution - 210 nm	<= 0.080	0.062
UV Absorbance of 0.25 M Solution - 200 nm	<= 0.200	0.116

For Laboratory, Research or Manufacturing Use  
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 17025:2005  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008, 17025:2005  
Panoli, India 9001:2008

Richard M Siberski  
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610