



Toluene

'BAKER ANALYZED'® HPLC Solvent
For Use in Liquid Chromatography

Product No. 9351
 Lot No. H41N54
 Release Date 11/16/2009

Certificate of Analysis

TEST	SPECIFICATION	RESULT
Assay (by GC, corrected for water)	99.7 % min.	99.9 %
Ultraviolet Absorbance		
(1.00-cm cell vs. water):		
350 nm	0.01 max.	0.005
300 nm	0.10 max.	0.10
288 nm	0.40 max.	0.34
UV Cut-off, nm	285 max.	284
Residue after Evaporation	3 ppm max.	0.3 ppm
Substances Darkened by H ₂ SO ₄	Passes Test	Passes Test
Water (by KF, coulometric)	0.02 % max.	0.005 %

For Laboratory, Research or Manufacturing Use

Country of Origin: USA

ISO

Phillipsburg, NJ 9001.2008, 14001.2004
 Paris, KY 9001.2008
 Mexico City, Mexico 9001.2008
 Deventer, The Netherlands 9001.2008, 14001.2004, 13485.2003
 Selangor, Malaysia 9001.2008
 Panoli, India 9001.2008
 Gliwice, Poland 9001.2008, 17025.2005

Marcy M. Matlosz

Marcy M. Matlosz
 Director of QA & Regulatory Affairs

For questions on this Certificate of Analysis please contact Technical Services at 1-800-582-2537 or 908-859-2151
 Mallinckrodt Baker, Inc. • 222 Red School Lane • Phillipsburg, NJ 08865 • Phone: 908.859.2151 • Fax: 908.859.6905