

Ether, Ultra Low Water
BakerDRY
Contains No Preservative



Material No.: 9250-12
Batch No.: 0000081438
Manufactured Date: 2014/06/02
Expiration Date: 2015/04/02

Certificate of Analysis

| Test | Specification | Result |
|---|---------------|---------|
| ACS – Assay ((C ₂ H ₅) ₂ O) (by GC) | >= 99.0 % | 100.0 |
| ACS – Alcohol (C ₂ H ₅ OH) | <= 0.05 % | < 0.01 |
| ACS – Carbonyl Compounds (as HCHO) (by polarography) | <= 0.001 % | < 0.001 |
| ACS – Color (APHA) | <= 10 | 5 |
| ACS – Peroxide (as H ₂ O ₂) | <= 1 ppm | < 1 |
| ACS – Residue after Evaporation | <= 5 ppm | < 1 |
| ACS – Titrable Acid (µeq/g) | <= 0.2 | 0.1 |
| ACS – Water (by KF, coulometric) | <= 10.0 ppm | 5.0 |

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Phillipsburg Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610