

p-(Dimethylamino)benzaldehyde
BAKER ANALYZED® A.C.S. Reagent

(Ehrlich's Reagent)



Material No.: J418-05
Batch No.: 0000107440
Manufactured Date: 2013/11/01
Retest Date: 2020/10/30

Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

Test	Specification	Result
Melting Point	73 – 75 °C	73
Solubility in HCl	Passes Test	PT
Solubility in Alcohol	Passes Test	PT
Color of HCl Solution	Passes Test	PT
Color of Alcohol Solution (APHA)	<= 60	<5
Residue after Ignition	<= 0.10 %	0.02

For Laboratory, Research or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610