

Hexafluoro-2-propanol  
BAKER ANALYZED® Reagent



Material No.: N151-05  
Batch No.: 0000110294  
Manufactured Date: 2014/03/15  
Retest Date: 2021/03/13

## Certificate of Analysis

Test	Specification	Result
Assay (C <sub>3</sub> H <sub>2</sub> F <sub>6</sub> O) (by GC, corrected for water)	>= 99 %	100
Color (APHA)	<= 10	5
Residue after Evaporation	<= 0.100 %	< 0.005
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	<= 0.05 %	0.01

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 17025:2005  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008, 17025:2005  
Panoli, India 9001:2008

  
Richard M Siberski  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610