

p-(Dimethylamino)benzaldehyde  
BAKER ANALYZED® A.C.S. Reagent

(Ehrlich's Reagent)



Material No.: J418-05  
Batch No.: 0000138889  
Manufactured Date: 2014/05/01  
Retest Date: 2021/04/29

## Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

Test	Specification	Result
Melting Point	73 – 75 °C	73
Solubility in HCl	Passes Test	PT
Solubility in Alcohol	Passes Test	PT
Color of HCl Solution	Passes Test	PT
Color of Alcohol Solution (APHA)	<= 60	<5
Residue after Ignition	<= 0.10 %	< 0.02

For Laboratory, Research or Manufacturing Use  
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008, 17025:2005  
Panoli, India 9001:2008

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610