

Hydrogen Peroxide, 30%
BAKER
(Stabilized)
Formerly C.P.



Material No.: 2189-01
Batch No.: 0000148979
Manufactured Date: 2016/07/19
Expiration Date: 2018/01/17

Certificate of Analysis

Test	Specification	Result
Assay (H ₂ O ₂)	29.0 - 32.0 %	31.6
Residue after Evaporation	<= 0.0300 %	< 0.0002
Nitrogen Compounds (as N)	<= 0.001 %	< 0.001
Phosphate (PO ₄)	<= 0.02 %	< 0.02
Sulfate (SO ₄)	<= 0.001 %	< 0.001
Chloride (Cl)	<= 5 ppm	< 5
ACS - Heavy Metals (as Pb)	<= 1 ppm	< 1
Trace Impurities - Iron (Fe)	<= 0.5 ppm	< 0.1

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008

Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610