

Hydrogen Peroxide, 30%  
BAKER  
(Stabilized)  
Formerly C.P.



Material No.: 2189-01  
Batch No.: 0000160725  
Manufactured Date: 2016/12/03  
Expiration Date: 2018/06/03

## Certificate of Analysis

Test	Specification	Result
Assay (H <sub>2</sub> O <sub>2</sub> )	29.0 - 32.0 %	31.6
Residue after Evaporation	<= 0.0300 %	< 0.0300
Nitrogen Compounds (as N)	<= 0.001 %	< 0.001
Phosphate (PO <sub>4</sub> )	<= 0.02 %	< 0.02
Sulfate (SO <sub>4</sub> )	<= 0.001 %	< 0.001
Chloride (Cl)	<= 5 ppm	< 5
ACS - Heavy Metals (as Pb)	<= 1 ppm	< 1
Trace Impurities - Iron (Fe)	<= 0.5 ppm	< 0.5

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 13485:2012  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008  
Panoli, India 9001:2008

Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610