

Tetrahydrofuran (Stabilized)
HPLC
For Use in Gel Permeation Chromatography
Contains Butylated Hydroxytoluene as a Preservative



Material No.: 9440-03
Batch No.: 0000159849
Manufactured Date: 2016/12/16
Expiration Date: 2019/12/13

Certificate of Analysis

| Test | Specification | Result |
|---|-----------------|--------|
| Assay (by GC, corrected for water, exclusive of preservative) | $\geq 99.5 \%$ | 100.0 |
| Color (APHA) | ≤ 10 | 5 |
| Peroxide (as H_2O_2) | $\leq 0.005 \%$ | 0.001 |
| Stabilizer (BHT) | 75 – 400 ppm | 93 |
| Water (by KF, coulometric) | $\leq 0.010 \%$ | 0.010 |

For Laboratory, Research or Manufacturing Use
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Country of Origin: US



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610