

Hydrogen Peroxide, 30%  
BAKER  
(Stabilized)  
Formerly C.P.



Material No.: 2189-01  
Batch No.: 0000195196  
Manufactured Date: 2018/03/02  
Expiration Date: 2019/08/31  
Revision No: 1

## Certificate of Analysis

Test	Specification	Result
Assay (H <sub>2</sub> O <sub>2</sub> )	29.0 - 32.0 %	31.7
Residue after Evaporation	<= 0.0300 %	0.0002
Nitrogen Compounds (as N)	<= 0.001 %	< 0.001
Phosphate (PO <sub>4</sub> )	<= 0.02 %	< 0.01
Sulfate (SO <sub>4</sub> )	<= 0.001 %	< 0.001
Chloride (Cl)	<= 5 ppm	< 1
ACS - Heavy Metals (as Pb)	<= 1 ppm	< 1
Trace Impurities - Iron (Fe)	<= 0.5 ppm	< 0.1

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003  
Gliwice, Poland 9001:2008, 13485:2012  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2008  
Panoli, India 9001:2008

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610