

Mannitol, Powder
BAKER ANALYZED® A.C.S. Reagent



Material No.: 2554-01
Batch No.: 0000199073
Manufactured Date: 2016/05/27
Retest Date: 2023/05/26
Revision No: 1

Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

Test	Specification	Result
ACS – Specific Rotation $[\alpha]^{25}_{\text{D}}$ (+)	23.3 – 24.3 Degree	23.8
ACS – Insoluble Matter	$\leq 0.01 \%$	< 0.01
ACS – Loss on Drying at 105°C	$\leq 0.05 \%$	0.02
ACS – Residue after Ignition	$\leq 0.01 \%$	< 0.01
ACS – Titrable Acid (meq/g)	≤ 0.0008	0.0001
ACS – Reducing Sugars	Passes Test	PT
Trace Impurities – ACS – Heavy Metals (as Pb)	$\leq 5 \text{ ppm}$	< 5

For Laboratory, Research or Manufacturing Use

Country of Origin: FR
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610