

Tetrahydrofuran (Stabilized)
HPLC
For Use in Gel Permeation Chromatography
Contains Butylated Hydroxytoluene as a Preservative



Material No.: 9440-03
Batch No.: 0000194395
Manufactured Date: 2018/03/09
Expiration Date: 2021/03/05
Revision No: 1

Certificate of Analysis

Test	Specification	Result
Assay (by GC, corrected for water, exclusive of preservative)	$\geq 99.5 \%$	99.9
Color (APHA)	≤ 10	5
Peroxide (as H_2O_2)	$\leq 0.005 \%$	< 0.005
Stabilizer (BHT)	75 – 400 ppm	99
Water (by KF, coulometric)	$\leq 0.010 \%$	< 0.001

For Laboratory, Research or Manufacturing Use
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway, Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610