

Methanol
ULTRA RESI-ANALYZED
For Purge and Trap Analysis



Material No.: 9077-02
Batch No.: 0000202813
Manufactured Date: 2018/05/21
Expiration Date: 2020/05/18
Revision No: 1

Certificate of Analysis

| Test | Specification | Result |
|---|---------------|--------|
| Assay (CH ₃ OH) (by GC, corrected for water) | >= 99.9 % | 100.0 |
| Residue after Evaporation | <= 1.0000 ppm | 0.2000 |
| Titration Acid (µeq/g) | <= 0.3 | 0.2 |
| Titration Base (µeq/g) | <= 0.1 | <0.01 |
| Water (by KF, coulometric) | <= 0.08 % | < 0.01 |
| Photoionization Detection (PID) Below CRQL | Passes Test | PT |
| Electroconductivity Detection (ELCD) Below CRQL | Passes Test | PT |

For Laboratory, Research or Manufacturing Use
Performance Tested for Use in EPA Methods
500 Series for Drinking Water
600 Series for Wastewater
846 for Solid Waste

Country of Origin: US
Packaging Site: Phillipsburg Mfg Ctr & DC



Phillipsburg, NJ 9001:2015, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Gliwice, Poland 9001:2008
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2015, 17025:2005
Panoli, India 9001:2015


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610