

Hexafluoro-2-propanol  
BAKER ANALYZED® Reagent



Material No.: N151-05  
Batch No.: 0000207077  
Manufactured Date: 2017/12/05  
Retest Date: 2024/12/03  
Revision No: 1

## Certificate of Analysis

Test	Specification	Result
Assay (C <sub>3</sub> H <sub>2</sub> F <sub>6</sub> O) (by GC, corrected for water)	>= 99 %	100
Color (APHA)	<= 10	5
Residue after Evaporation	<= 0.100 %	0.001
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	<= 0.05 %	0.01

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2015, FSSC 22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Gliwice, Poland 9001:2008  
Selangor, Malaysia 9001:2008  
Dehradun, India 9001:2015, 14001:2015, 13485:2016  
Mumbai, India 9001:2015  
Panoli, India 9001:2015

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600  
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610