

Ether, Ultra Low Water  
BakerDRY  
Contains No Preservative



Material No.: 9250-12  
Batch No.: 0000213261  
Manufactured Date: 2018/06/18  
Expiration Date: 2019/04/18  
Revision No: 1

## Certificate of Analysis

Test	Specification	Result
ACS - Assay ((C <sub>2</sub> H <sub>5</sub> ) <sub>2</sub> O) (by GC)	>= 99.0 %	99.9
ACS - Alcohol (C <sub>2</sub> H <sub>5</sub> OH)	<= 0.05 %	< 0.01
ACS - Carbonyl Compounds (as HCHO) (by polarography)	<= 0.001 %	< 0.001
ACS - Color (APHA)	<= 10	<5
ACS - Peroxide (as H <sub>2</sub> O <sub>2</sub> )	<= 1 ppm	< 1
ACS - Residue after Evaporation	<= 5 ppm	< 2
ACS - Titrable Acid (µeq/g)	<= 0.2	< 0.1
ACS - Water (by KF, coulometric)	<= 10.0 ppm	7.1

For Laboratory, Research or Manufacturing Use

Country of Origin: US  
Packaging Site: Phillipsburg Mfg Ctr & DC



Phillipsburg, NJ 9001:2015, FSC22000  
Paris, KY 9001:2008  
Mexico City, Mexico 9001:2008  
Gliwice, Poland 9001:2015, 13485:2012  
Selangor, Malaysia 9001:2008  
Dehradun, India, 9001:2008, 14001:2004, 13485:2003  
Mumbai, India, 9001:2015, 17025:2005  
Panoli, India 9001:2015

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700