

Tetrahydrofuran (Stabilized)
HPLC
For Use in Gel Permeation Chromatography
Contains Butylated Hydroxytoluene as a Preservative



Material No.: 9440-03
Batch No.: 0000214351
Manufactured Date: 2018/10/22
Expiration Date: 2021/10/18
Revision No: 1

Certificate of Analysis

Test	Specification	Result
Assay (by GC, corrected for water, exclusive of preservative)	$\geq 99.5 \%$	100.0
Color (APHA)	≤ 10	< 5
Peroxide (as H_2O_2)	$\leq 0.005 \%$	0.003
Stabilizer (BHT)	75 – 400 ppm	214
Water (by KF, coulometric)	$\leq 0.010 \%$	0.005

For Laboratory, Research or Manufacturing Use
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2015, FSC22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Gliwice, Poland 9001:2015, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2015, 17025:2005
Panoli, India 9001:2015


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700