

p-(Dimethylamino)benzaldehyde  
BAKER ANALYZED® A.C.S. Reagent

(Ehrlich's Reagent)



Material No.: J418-05  
Batch No.: 0000226880  
Manufactured Date: 2017/06/01  
Retest Date: 2024/05/30  
Revision No: 1

## Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

Test	Specification	Result
Melting Point	73 – 75 °C	73
Solubility in HCl	Passes Test	PT
Solubility in Alcohol	Passes Test	PT
Color of HCl Solution	Passes Test	PT
Color of Alcohol Solution(APHA)	<= 60	<5
Residue after Ignition	<= 0.10 %	< 0.10

For Laboratory, Research or Manufacturing Use  
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB  
Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700