

Tetrahydrofuran (Stabilized)
HPLC
For Use in Gel Permeation Chromatography
Contains Butylated Hydroxytoluene as a Preservative



Material No.: 9440-03
Batch No.: 0000068669
Manufactured Date: 2014/03/07
Expiration Date: 2017/03/03

Certificate of Analysis

| Test | Specification | Result |
|---|-----------------|-----------|
| Assay (by GC, corrected for water, exclusive of preservative) | $\geq 99.5 \%$ | 100.0 |
| Color (APHA) | ≤ 10 | 5 |
| Peroxide (as H_2O_2) | $\leq 0.005 \%$ | < 0.001 |
| Stabilizer (BHT) | 75 – 400 ppm | 200 |
| Water (by KF, coulometric) | $\leq 0.010 \%$ | 0.007 |

For Laboratory, Research or Manufacturing Use
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Country of Origin: US



Phillipsburg, NJ 9001:2008, 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610