

Barium Sulfate, Powder
BAKER ANALYZED® Reagent



Material No.: 1030-05
Batch No.: 0000110791
Manufactured Date: 2014/07/26
Retest Date: 2021/07/24

Certificate of Analysis

Test	Specification	Result
Assay (BaSO ₄)	>= 97.0 %	99.7
Acidity (as H ₂ SO ₄)	<= 0.01 %	< 0.01
Water-Soluble Salts	<= 0.02 %	< 0.02
Chloride (Cl)	<= 0.005 %	0.005
Phosphate (PO ₄)	<= 0.001 %	< 0.001
Heavy Metals (as Pb)	<= 0.001 %	< 0.001
Trace Impurities – Iron (Fe)	<= 0.001 %	< 0.001
Trace Impurities – Arsenic (As)	<= 1.000 ppm	< 1.000

For Laboratory, Research or Manufacturing Use

Country of Origin: DE
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610