

1-Octanesulfonic Acid, Sodium Salt
HPLC
For Ion-Pair Chromatography of Basic Compounds



Material No.: 2818-05
Batch No.: 0000166830
Manufactured Date: 2016/12/01
Retest Date: 2023/11/30

Certificate of Analysis

| Test | Specification | Result |
|--|---------------|--------|
| Assay (as C ₈ H ₁₇ SO ₃ Na) (by acidimetry) | >= 98.0 % | 100.1 |
| UV Absorbance of 0.25 M Solution – 250 nm | <= 0.050 | 0.010 |
| UV Absorbance of 0.25 M Solution – 240 nm | <= 0.050 | 0.010 |
| UV Absorbance of 0.25 M Solution – 230 nm | <= 0.050 | 0.014 |
| UV Absorbance of 0.25 M Solution – 220 nm | <= 0.060 | 0.024 |
| UV Absorbance of 0.25 M Solution – 210 nm | <= 0.080 | 0.058 |
| UV Absorbance of 0.25 M Solution – 200 nm | <= 0.200 | 0.085 |

For Laboratory, Research or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004, FSSC 22000
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 13485:2012
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008
Panoli, India 9001:2008

Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor Performance Materials, LLC.

3477 Corporate Parkway. Center Valley, PA 18034. U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610