

Ether, Ultra Low Water
BakerDRY
Contains No Preservative



Material No.: 9250-12
Batch No.: 0000044475
Manufactured Date: 2013/05/03
Expiration Date: 2014/05/03

Certificate of Analysis

Test	Specification	Result
ACS - Assay ((C ₂ H ₅) ₂ O) (by GC)	>= 99.0 %	99.9
ACS - Alcohol (C ₂ H ₅ OH)	<= 0.05 %	< 0.01
ACS - Carbonyl Compounds (as HCHO) (by polarography)	<= 0.001 %	< 0.001
ACS - Color (APHA)	<= 10	<5
ACS - Peroxide (as H ₂ O ₂)	<= 1 ppm	< 1
ACS - Residue after Evaporation	<= 5 ppm	< 2
ACS - Titrable Acid (µeq/g)	<= 0.2	< 0.1
ACS - Water (by KF, coulometric)	<= 10.0 ppm	9.0

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Phillipsburg Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610