

p-(Dimethylamino)benzaldehyde
BAKER ANALYZED® A.C.S. Reagent

(Ehrlich's Reagent)



Material No.: J418-05
Batch No.: 0000039965
Manufactured Date: 2011/02/01
Retest Date: 2016/01/31

Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

| Test | Specification | Result |
|----------------------------------|---------------|--------|
| Melting Point | 73 – 75 °C | 73 |
| Solubility in HCl | Passes Test | PT |
| Solubility in Alcohol | Passes Test | PT |
| Color of HCl Solution | Passes Test | PT |
| Color of Alcohol Solution (APHA) | <= 60 | <5 |
| Residue after Ignition | <= 0.10 % | 0.03 |

For Laboratory, Research or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB
Packaging Site: Paris Mfg Ctr & DC



Phillipsburg, NJ 9001:2008, 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, The Netherlands 9001:2008, 14001:2004, 13485:2003
Gliwice, Poland 9001:2008, 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India, 9001:2008, 14001:2004, 13485:2003
Mumbai, India, 9001:2008, 17025:2005
Panoli, India 9001:2008

Richard M Siberski
Global Director of Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.573.2600
Avantor™ Performance Materials Inc.

3477 Corporate Parkway, Suite #200, Center Valley, PA 18034, U.S.A. Phone: 610.573.2600 . Fax: 610.573.2610