

1-Butanol
PHOTREX® Reagent
For Spectrophotometry
(butyl alcohol)



Material No.: 9189-01
Batch No.: 0000235376
Manufactured Date: 2019/06/15
Expiration Date: 2021/06/12
Revision No: 1

Certificate of Analysis

Test	Specification	Result
ACS – Assay ($\text{CH}_3\text{CH}_2\text{CH}_2\text{CH}_2\text{OH}$) (by GC, corrected for water)	$\geq 99.4 \%$	100.0
ACS – Color (APHA)	≤ 10	< 5
ACS – Residue after Evaporation	$\leq 0.0005 \%$	0.0004
ACS – Titrable Acid ($\mu\text{eq/g}$)	≤ 0.8	0.6
ACS – Butyl Ether ($\text{C}_8\text{H}_{18}\text{O}$)(by GC)	$\leq 0.2 \%$	< 0.1
Carbonyl Compounds – Butyraldehyde	$\leq 0.01 \%$	< 0.01
Water (by KF, coulometric)	$\leq 0.05 \%$	0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) – 400–300 nm	≤ 0.01	< 0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) – 280 nm	≤ 0.02	0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) – 260 nm	≤ 0.04	0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) – 240 nm	≤ 0.10	0.04
Ultraviolet Absorbance (1.00-cm cell vs. water) – 220 nm	≤ 0.50	0.23
Ultraviolet Absorbance (1.00-cm cell vs. water) – 215 nm	≤ 1.00	0.34

For Laboratory, Research or Manufacturing Use

Product Information (not specifications):

Boiling Point (typical)..... 117.7 °C.

Windows of Infrared Transmittance (0.1-mm path, 50–100% T), μm :

2.0 – 2.8 ; 4.3 – 6.0

Material No.: 9189-01
Batch No.: 0000235376

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700