

Methanol
ULTRA RESI-ANALYZED
For Purge and Trap Analysis



Material No.: 9077-02
Batch No.: 0000230446
Manufactured Date: 2019/04/11
Expiration Date: 2022/04/10
Revision No: 1

Certificate of Analysis

| Test | Specification | Result |
|---|---------------|--------|
| Assay (CH ₃ OH) (by GC, corrected for water) | >= 99.9 % | 100.0 |
| Residue after Evaporation | <= 1.0000 ppm | 0.3000 |
| Titration Acid (µeq/g) | <= 0.3 | 0.2 |
| Titration Base (µeq/g) | <= 0.1 | <0.01 |
| Water (by KF, coulometric) | <= 0.08 % | < 0.01 |
| Photoionization Detection (PID) Below CRQL | Passes Test | PT |
| Electroconductivity Detection (ELCD) Below CRQL | Passes Test | PT |

For Laboratory, Research or Manufacturing Use
Performance Tested for Use in EPA Methods
500 Series for Drinking Water
600 Series for Wastewater
846 for Solid Waste

Country of Origin: US
Packaging Site: Phillipsburg Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700