

Tetrahydrofuran (Stabilized)  
HPLC  
For Use in Gel Permeation Chromatography  
Contains Butylated Hydroxytoluene as a Preservative



Material No.: 9440-03  
Batch No.: 0000249999  
Manufactured Date: 2020/01/23  
Expiration Date: 2023/01/22  
Revision No: 1

## Certificate of Analysis

Test	Specification	Result
Assay (by GC, corrected for water, exclusive of preservative)	$\geq 99.5 \%$	100.0
Color (APHA)	$\leq 10$	$< 5$
Peroxide (as $\text{H}_2\text{O}_2$ )	$\leq 0.005 \%$	$< 0.001$
Stabilizer (BHT)	75 – 400 ppm	258
Water (by KF, coulometric)	$\leq 0.010 \%$	0.004

For Laboratory, Research or Manufacturing Use  
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Country of Origin: US  
Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700