

p-(Dimethylamino)benzaldehyde  
BAKER ANALYZED® A.C.S. Reagent

(Ehrlich's Reagent)



Material No.: J418-05  
Batch No.: 0000250637  
Manufactured Date: 2019/11/01  
Retest Date: 2026/10/30  
Revision No: 1

## Certificate of Analysis

Meets ACS Reagent Chemical Requirements,

| Test                             | Specification | Result |
|----------------------------------|---------------|--------|
| Melting Point                    | 73 – 75 °C    | 73     |
| Solubility in HCl                | Passes Test   | PT     |
| Solubility in Alcohol            | Passes Test   | PT     |
| Color of HCl Solution            | Passes Test   | PT     |
| Color of Alcohol Solution (APHA) | <= 60         | < 60   |
| Residue after Ignition           | <= 0.10 %     | < 0.10 |

For Laboratory, Research or Manufacturing Use  
Meets Reagent Specifications for testing USP/NF monographs

Country of Origin: GB  
Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700