

Tetrahydrofuran (Stabilized)
HPLC
For Use in Gel Permeation Chromatography
Contains Butylated Hydroxytoluene as a Preservative



Material No.: 9440-03
Batch No.: 0000203777
Manufactured Date: 2018/08/07
Expiration Date: 2021/08/03
Revision No: 2

Certificate of Analysis

Test	Specification	Result
Assay (by GC, corrected for water, exclusive of preservative)	$\geq 99.5 \%$	100.0
Color (APHA)	≤ 10	< 5
Peroxide (as H ₂ O ₂)	$\leq 0.005 \%$	0.004
Stabilizer (BHT)	75 – 400 ppm	255
Water (by KF, coulometric)	$\leq 0.010 \%$	0.009

For Laboratory, Research or Manufacturing Use
Stabilized with approximately 0.01% Butylated Hydroxytoluene

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700