

Tetrahydrofuran
HPLC
For use in Liquid Chromatography (HPLC & UHPLC) &
Spectrophotometry
Contains No Preservative



Material No.: 9441-33
Batch No.: 0000274892
Manufactured Date: 2021/02/04
Expiration Date: 2022/02/03
Revision No: 1

Certificate of Analysis

Test	Specification	Result
Ultraviolet Absorbance (1.00-cm cell vs. water) – 400–320 nm	<= 0.01	< 0.01
Ultraviolet Absorbance (1.00-cm cell vs. water) – 280 nm	<= 0.030	0.112
Ultraviolet Absorbance (1.00-cm cell vs. water) – 254 nm	<= 0.15	0.14
Ultraviolet Absorbance (1.00-cm cell vs. water) – 230 nm	<= 0.45	1.16
Ultraviolet Absorbance (1.00-cm cell vs. water) – UV Cut-off, nm	<= 212	232
Assay (by GC, corrected for water)	>= 99.8 %	100.0
Peroxide (as H ₂ O ₂) (atTime of Fill)	<= 0.005 %	0.002
Residue after Evaporation	<= 2.0 ppm	124.5
Water (by KF, coulometric)	<= 0.02 %	0.01

For Laboratory, Research or Manufacturing Use

Country of Origin: US
Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone: 610.386.1700